



# ***Brava® Protective Ring Convex***

Customised convexity for increased protection

***Brava® Protective Ring Convex***

# The first combined solution designed for the outward peristomal body profile

How many people can benefit from this solution?

38 %<sup>3</sup>

X

29 %<sup>3</sup>

=

11%

of Canadians with ostomies have an outward peristomal body profile.

of those have a flush or retracted stoma.

of Canadians with ostomies.

Together, Brava<sup>®</sup> Protective Ring Convex and SenSura<sup>®</sup> Mio Concave achieved 99% skin contact area\*

**1 SenSura Mio Concave**


The only ostomy product designed for bulges and curves, SenSura Mio Concave is proven to provide better skin contact and reduce risk of leakage compared to flat baseplates.

**2 Brava Protective Ring Convex**

At the same time, Brava Protective Ring Convex is designed to provide convexity to the peristomal area, with double adhesive layers for added protection against leakage.

**3 Double protection for flush and retracted stomas**

Together, Brava Protective Ring Convex and SenSura Mio Concave provide a snug and secure seal between the baseplate and the user's skin.



\*From bench testing simulating outward peristomal body profiles with flush or retracted stomas. Actual results may vary based on individual situation.

# Some stomas need a little extra support

For anyone who lives with a stoma, leakage from their ostomy product can be a major cause of anxiety – as well as causing skin complications.

Brava® Protective Ring is designed to provide a better fit between the skin and the baseplate and is proven to significantly reduce leakage compared to other rings on the market<sup>1</sup>.

However, some stomas need a little extra support and that's why we created **Brava® Protective Ring Convex**. It provides customised convexity for patients with flush and retracted stomas – helping them to live active, healthy lives.



## 81%

of users prefer Brava® Protective Ring over their current seal.<sup>2</sup>



## New Brava® Protective Ring Convex

The new Brava Protective Ring Convex is designed to protect the skin against leakage for flush or retracted stomas and patients with folds, creases or scarring around their stoma.

### Brava® Protective Ring adhesive

- High tack and fills small creases
- Skin friendly adhesive

### SenSura® Mio adhesive


- Erosion resistant adhesive towards the stoma
- High flexibility for protective seal during body movement

### Designed for a perfect fit

- Steep sides designed to provide convexity to peristomal area
- Thin edges fit snugly around peristomal area

### Shapable with sustained hold

- Shapable to fit individual stoma shapes, including oval and loop
- Sustains convexity after shaping



# Customised convexity for increased protection

Brava® Protective Ring Convex provides customised convexity for people with a flush or retracted stoma. The steep sides are designed to provide convexity to the peristomal area, while the double adhesive layers help protect against leakage and at the same time, is designed to reduce the moisture of the skin to maintain a good pH balance.



## How the Brava Protective Ring Convex works

1


Brava Protective Ring Convex is placed around the stoma.

2

The double adhesive layers adhere securely to the skin.

3

The shape is designed to provide convexity to the peristomal area.




# Customise your solution

Brava® Protective Ring Convex comes in an extensive range of sizes that can be matched to individual needs. It can also be formed around all stoma shapes, while sustaining it's convexity.


## Round or oval stoma?

Shapable for different stomas, including into ovals for loop or temporary stomas.




## Mix and match

5 mm intervals provide an easy to match sizing with SenSura Mio Concave. Can also be used in combination with flat baseplates.



## Sizing for the perfect fit

Brava Protective Ring Convex is available in six sizes with 5 mm intervals.



	20 mm	25 mm	30 mm	35 mm	40 mm	45 mm
Brava Protective Ring Convex	12090	12091	12092	12093	12094	12095

1. BPS Leakage Test: Coloplast Lab Report. 2020 (tested against Dansac TRE Seals, Eakin Cohesive Seals, Hollister CeraRings, Salts Mouldable Seals with Aloe, Welland Hyperseals® with Manuka Honey, ConvaTec Stomahesive® Seals). Data on file. VV-0283280.
2. Coloplast, Product Evaluation, Brava Protective Seal, 2017, N=281 users, Data-on-File (VV-0204884)
3. Coloplast, Market Research, Ostomy Life Study (OLS) 2019. Data on File.

[Ostomy Care](#) / [Continence Care](#) / [Wound & Skin Care](#) / [Interventional Urology](#)

Coloplast Canada, 2401 Bristol Circle, Suite A205, Oakville, ON L6H 5S9

[www.coloplast.ca](http://www.coloplast.ca) The Coloplast logo is a registered trademark of Coloplast.

© 2021 Coloplast Canada. All rights reserved.

